

# Easy to use



*1. Clean the site of bleeding*



*2. Fill the site of bleeding with light pressure either with intact sponge or cut it to the desired size*



*3. Allow for 2-5 minutes for haemostasis*

Manufactured by  
**Goodwill Lifesciences,**

157-158, RJD Integrated Textile Park,  
Paal - Hazira Road, Ichhapore, Surat - 394 510, Gujarat.

[www.goodwill-lifesciences.com](http://www.goodwill-lifesciences.com)

e-mail: [info@goodwill-lifesciences.com](mailto:info@goodwill-lifesciences.com),

[hemosponge@gmail.com](mailto:hemosponge@gmail.com)

Tel.: 0261-2974173, M: 09879744173

**Mfg Lic. No.: G/25/2165**

**Goodwill**  
Patient Safety First

## HEMOSPONGE

Dental cubes SL/XL

Dental  
Haemostatic  
Sponge



**Fast**

Haemostasis in 2-4 minutes

**Safe**

Non-toxic, Non-allergic and biocompatible

**Effective**

Stop all kinds of oozing blood in dental and oral surgeries

**Highly absorbent**

Absorbs 35 – 40 times of its own weight of blood/fluids

**pH neutral**

Can be soaked in variety of pharmaceuticals/biologics

**Biodegradable**

Manufactured from pharma gelatin

**Absorbable**

Absorb in body in 2-3 days

**Gamma sterile**

No need to sterilize before use

**Quality assurance**

Complies US Pharmacopoeia standards

**Available in two sizes**

Dental SL: 10 × 10 × 10 mm

Dental XL: 20 × 20 × 07 mm

**Effective haemostatic for**

- 1) Dental surgeries
  - Endodontic
  - Periodontics
  - Exodontia
  - Prosthodontics
- 2) Tissue biopsy
- 3) Oral and maxillofacial surgery

**Mechanism of action**

HEMOSPONGE arrests the bleeding by the formation of an artificial clot. In addition to serving as a mechanical obstruction to bleeding surfaces, the porous structure of HEMOSPONGE facilitates the process of blood clotting. At the moment thrombocytes (platelets) come in contact with the 3D matrix of the sponge, they change their surface characters to release a thromboplastin, activating the cascade of clotting.